

MOLECULAR BIOLOGY 2012 **CATALOG**

GemTaqTM DNA Polymerase

GemTaqTM DNA Polymerase is a special formulation of standard Taq DNA Polymerase that features significantly improved amplifications of DNA fragments. GemTaqTM DNA Polymerase has an elongation rate higher than regular Taq DNA Polymerase, resulting in robust yields and greater sensitivity of PCR amplifications within a short reaction time.

Advantages


- Higher efficiency and greater sensitivity in a wide range of PCR assays when compared to a standard Taq DNA Polymerase.
- Suited for amplifications of long DNA fragments.
- Shorter reaction times due to higher elongation rate.
- Unbeatable value in terms of cost per unit and performance.

Human genomic DNA (ng/25 μl) 12 6 3 1.5 0.75 0.37 0.19 0.1 GemTaq Taq 1 Taq 2

The $\beta\text{-Actin}$ gene fragment (850 bp) was amplified in 25 μl from the indicated amounts of human genomic DNA as a template. Each PCR reaction was performed using 0.25 U of GemTaq^TM DNA Polymerase and standard Taq DNA Polymerases (Taq 1 from competitor 1, Taq 2 from competitor 2) in supplied Reaction Buffers, respectively during 33 cycles.

Applications

- Routine and high-throughput PCR.
- Primer extension.
- TA cloning.





Lambda DNA fragments (3000 bp and 5000 bp) were amplified in 50 μ l with 0.5 U of GemTaqTM DNA Polymerase and standard Taq DNA polymerases in supplied Reaction Buffers, respectively. Each PCR reaction was performed from 5 ng of Lambda DNA as a template during 28 cycles.

Product	Catalog #	Size	Conc.	Price USD
GemTaq TM DNA Polymerase and 5xReaction Buffer (with 7.5 mM Mg ²⁺)	EP012 EP013 EP014 EP015	400 units 2,000 units 5,000 units 15,000 units	5 u/μl 5 u/μl 5 u/μl 5 u/μl	\$ 45.00 \$ 175.00 \$ 400.00 \$ 1,050.00
GemTaq TM DNA Polymerase and Colored 5xReaction Buffer (with 7.5 mM Mg ²⁺)	EP022 EP023 EP024 EP025	400 units 2,000 units 5,000 units 15,000 units	5 u/μl 5 u/μl 5 u/μl 5 u/μl	\$ 45.00 \$ 175.00 \$ 400.00 \$ 1,050.00
GemTaq [™] 5xMaster Mix	EPM11 EPM12 EPM13 EPM14	0.5 ml (100 reactions) 2.5 ml (500 reactions) 5.0 ml (1,000 reactions) 10.0 ml (2,000 reactions)		\$ 25.00 \$ 100.00 \$ 180.00 \$ 320.00
GemTaq [™] Colored 5xMaster Mix	EPM21 EPM22 EPM23 EPM24	0.5 ml (100 reactions) 2.5 ml (500 reactions) 5.0 ml (1,000 reactions) 10.0 ml (2,000 reactions)		\$ 25.00 \$ 100.00 \$ 180.00 \$ 320.00

Colored 5xReaction Buffer

The Colored 5xReaction Buffer for amplifications was specially formulated to provide direct loading of the PCR samples onto gels without the need to add a loading sample buffer. It contains a compound that increases sample density for electrophoresis and two dyes (Red and Yellow) that separate during electrophoresis to allow for easy monitoring. The Red and Yellow dyes migrate at the same rate as 800-1000 bp DNA fragments and 20-30 bp DNA fragments in a 1% agarose gel, respectively.





Lambda DNA fragment (5000 bp) was amplified in 50 μ l with 2U (lane 1), 1.25U (lane 2), 0.83U (lane 3), 0.55U (lane 4), 0.37U (lane 5) and 0.25U (lane 6) of GemTaqTM DNA Polymerase in the Colored Reaction Buffer. Each PCR reaction was performed from 2 ng of Lambda DNA as a template during 30 cycles. Lane M: 1 kb DNA ladder.

GemTaq™ Hot Start DNA Polymerase

GemTaqTM Hot Start DNA Polymerase is a heat-activated high-performance proprietary formulation of GemTaqTM DNA Polymerase. The modified polymerase activity is quickly restored during the initial denaturation step allowing for a hot-start PCR. Therefore GemTaqTM Hot Start DNA Polymerase provides improved specificity and very high PCR sensitivity by eliminating the presence of non-specifics products.

Advantages

- · Ultra-high specificity.
- · Fast activation.
- Recommended for low copy assays.
- Amplifies fragments up to 5 kb.
- Cost-effective.

Applications

- Hot start PCR.
- Low-copy target.
- · High throughput PCR applications.
- Complex genomic templates.
- · TA cloning.

Product	Catalog #	Size	Conc.	Price USD
GemTaq TM Hot Start DNA Polymerase and $5xReaction$ Buffer (with 7.5 mM Mg^{2+})	EP032	400 units	5 u/μl	\$ 55.00
	EP033	2,000 units	5 u/μl	\$ 220.00
	EP034	5,000 units	5 u/μl	\$ 500.00
	EP035	15,000 units	5 u/μl	\$ 1,300.00
GemTaq™ Hot Start DNA Polymerase and Colored 5xReaction Buffer (with 7.5 mM Mg ²⁺)	EP042	400 units	5 u/μl	\$ 55.00
	EP043	2,000 units	5 u/μl	\$ 220.00
	EP044	5,000 units	5 u/μl	\$ 500.00
	EP045	15,000 units	5 u/μl	\$ 1,300.00
GemTaq [™] Hot Start 5xMaster Mix	EPM31	0.5 ml (100 reactions)		\$ 30.00
	EPM32	2.5 ml (500 reactions)		\$ 130.00
	EPM33	5.0 ml (1,000 reactions)		\$ 225.00
	EPM34	10.0 ml (2,000 reactions)		\$ 400.00
GemTaq™ Hot Start Colored 5xMaster Mix	EPM41 EPM42 EPM43 EPM44	0.5 ml (100 reactions) 2.5 ml (500 reactions) 5.0 ml (1,000 reactions) 10.0 ml (2,000 reactions)		\$ 30.00 \$ 130.00 \$ 225.00 \$ 400.00

SNPase Hot Start DNA Polymerase

SNPase Hot Start is a thermostable DNA polymerase designed for a single-nucleotide polymorphism (SNP) detection. This enzyme has only 5'-3' polymerase activity and lacks exonuclease activity. SNPase Hot Start DNA Polymerase provides very high enzymatic specificity and high fidelity incorporation of dNTPs and ddNTPs that is 10-15 times higher than a regular *Taq* DNA Polymerase. The enzyme has a modified polymerase activity that is restored during the initial denaturation step allowing for a hot-start PCR and eliminates the presence of non-specifics products. **SNPase was recognized as one of the best Polymerases in SNP detection by allele-specific PCR and mini-sequencing.**

Advantages


- High sensitivity to mismatches at 3'-end of primers.
- High fidelity of dNTPs and ddNTPs incorporation.
- Lack of 5'-3' and 3'-5' exonuclease activity.
- High specificity, low level of background PCR.
- Fast activation (5-10 sec, 95°C).
- TA cloning.

Applications

- SNP- genotyping, allele specific PCR, allele specific primer extension, mini-sequencing.
- High fidelity PCR (up to 1000 bp).
- Multiplex PCR.
- Real Time PCR with intercalating dyes (SYBR Green, Eva Green).

Product	Catalog #	Size	Conc.	Price USD
SNPase Hot Start DNA Polymerase and 5xReaction Buffer (with 12.5 mM Mg ²⁺)	EP052	400 units	12 u/μl	\$ 60.00
	EP053	2,000 units	12 u/μl	\$ 240.00
	EP054	5,000 units	12 u/μl	\$ 520.00





MGQuest is a fast growing private American company located in the beautiful state of Washington.

The company offers products and services for research in life sciences. Our goal is to provide research communities with the highest quality products at very competitive prices.

We support advances in science by inventing and marketing products that accelerate and improve research.



web www.mgquest.com

Email sales@mgquest.com Email techserv@mgquest.com

> USA P.O. Box 5739 Lynnwood WA 98046-5739